

# SpotOptics

The software people for optics

## OPTINO

### VERSATILE WAVEFRONT SENSOR

- Accurate metrology in single and double pass
- Lenses, mirrors and laser beams
- Any focal length and diameter
- Large dynamic range
- Adaptable for production
- Wide wavelength range with different cameras (193nm-10.6 $\mu$ )
- Different lenslets available according to application



## Optino: models



### OPTINO-MU

Small Footprint. Max: 45x45 spots  
 USB3  
 Wavelength range: 380nm-1100nm  
 Readout speed: 90Hz  
 Processing speed: 15Hz  
 37(H) x 30(W) x 79(L) mm



### OPTINO-UV-VIS

Max 75x75 spots  
 Thinned CCD chip  
 Gigabit Ethernet  
 Wavelength range: 193nm-1100nm  
 Readout speed: 7.5Hz  
 Processing speed: 7.5Hz





### OPTINO-NIR

Max 45x45 spots  
 InGaAs camera  
 Gigabit Ethernet  
 Wavelength range: 950nm-1700nm  
 Readout speed: 100Hz  
 Processing speed: 50Hz

## Technical Specifications: Hardware

<b>Focal ratios covered</b>	Practically any focal ratio can be tested with appropriate collimator
<b>Diameter of optical element that can be tested in parallel light</b>	Practically any diameter can be tested with appropriate accessories
<b>Test setup</b>	Double or single pass
<b>Diameter and focal length of standard lenslets</b>	(0.2mm, 11mm), (0.2mm, 22mm), (0.3mm, 41mm). Others: (0.13mm, 7mm). Special
<b>Camera - connections</b>	Cameras with Gigabit and USB2 and USB3 connection available. 10-bits to 16-bits
<b>No of spots (0.2mm diameter lenslet array and USB3 camera). 9mm pupil</b>	45x45 with (with 0.2mm array). 69x69 spots (with 0.13mm array)
<b>Maximum no. of spots (array with Gigabit camera). 15mm pupil</b>	75x75 (with 0.2mm array). 115x115 spots (with 0.13mm array)
<b>Measure aspherical elements</b>	~15% variation in longitudinal spherical aberration
<b>RMS repeatability of computation of Zernike polynomials</b>	1-2nm rms ( $\lambda/600$ - $\lambda/300$ ) at 633nm
<b>RMS repeatability of wavefront measurements</b>	$< \lambda/200$
<b>Accuracy</b>	$\lambda/10$ - $\lambda/100$ – depending on the accuracy of the calibration elements
<b>Dynamic range of measurements sub-pupil (tilt subtracted)</b>	$\pm 50 \lambda$
<b>Wavelength range</b>	UV (0.193-1.1 $\mu$ ), Vis (0.193-1.1 $\mu$ ), 0.95-1.7 $\mu$
<b>Light sources</b>	LEDs, LDs and Halogen lamps available (at one wavelength or in white light)
<b>Software (control and analysis)</b>	Sensoft
<b>Acquisition speed</b>	15-500Hz (camera dependent)
<b>Processing speed</b>	5-150 Hz (camera and PC dependent)
<b>Power requirement for collimator stepper motor (motorized version)</b>	12V/2A DC
<b>Trigger</b>	TTL 5V

## SENSOFT: THE SOFTWARE




### Sensoft: The modular software package

- ✓ Hardware control of Optino (camera, light source and motor)
- ✓ Full Shack-Hartmann (SH) analysis
- ✓ Zernike aberration coefficients: Annular, Standard and Fringe
- ✓ Aberration coefficients: Seidel
- ✓ Wavefront: Zonal and modal
- ✓ MTF, PSF and EE
- ✓ Strehl ratio
- ✓ Spot diagram and
- ✓ Diagnostics (alignment and best focal plane)
- ✓ Loop mode for on-line alignment and best focal plane) using coma, astigmatism and spherical aberration
- ✓ Loop mode for on-line adjustment of spacing of elements to get best focal plane using spherical aberration
- ✓ Average of SH images to decrease noise
- ✓ Automatic control of exposure time

### Optino in your production line:


- ✓ Optino – with its own PC - can easily be adapted to the production line
- ✓ Work in a closed-loop with the PC of the manufacturing machine
- ✓ A software module defines the communication protocol and transfers the results between the control PC of Optino and the production PC, minimizing cost and time to develop special software to install on user PC

# ON-LINE ALIGNMENT IN A FAST LOOP




- The alignment of complex optical systems becomes easy by monitoring coma and astigmatism in a continuous loop
- The individual (x, y) components of coma and astigmatism, as well as the total coefficients are displayed
- The optimization can be done for one component at a time, as the software can display one component of interest
- Optimal alignment is reached when the coma and astigmatism

# MTF MEASUREMENTS



MTF after subtracting the contributions of tilt and defocus present in the data.

## LAYOUT



## PHYSICAL

Dimensions: 17 (L) x 10 (W) x 10 (H) cm

Weight: 3 Kg

Cameras: USB3, Gigabit Ethernet

Motor power supply: 12V, 2A

## KEY FEATURES

### Measurement technique

Shack-Hartmann wavefront sensor

**Test in parallel light or at the lens focus**

Parallel light (with a collimator)

At the focus of the lens (with a pinhole)

Light sources with different wavelength available

### Calibration

High-quality parallel light source

Pinhole calibration unit

$\lambda/20$  Spherical and flat reference mirrors

### Versatile

Test any element in double or single pass

### Accessories

Full set of accessories available

## SOFTWARE

Easy alignment of lens group via software: the software gives graphical indication of the misalignment of the optical system, using coma and astigmatism.